

NEDD4 Antibody



PACO10807

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NEDD4

Uniprot

P46934

Synonyms:

neural precursor cell expressed, developmentally down-regulated 4; NEDD4; KIAA0093; MGC176705; NEDD4-1; RPF1 ; NEDD-4;

Immunogen:

Human NEDD4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A