

NEFH Antibody



PACO10810

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NEFH

Uniprot

P12036

Synonyms:

neurofilament, heavy polypeptide; NEFH; NFH ;

Immunogen:

Human NEFH.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A