

NFATC2IP Antibody



PACO10827

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NFATC2IP

Uniprot

Q8NCF5

Synonyms:

nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 2 interacting protein; NFATC2IP; FLJ14639; MGC126790; MGC138387 ; NIP45

Immunogen:

Human NFATC2IP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A