

NFIB Antibody



PACO10832

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NFIB

Uniprot

O00712

Synonyms:

nuclear factor I/B; NFIB; HMGIC/NFIB; NFI-RED; NFIB2; NFIB3 ;

Immunogen:

Human NFIB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A