

NFKBIA Antibody



PACO10836

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NFKBIA

Uniprot

P25963

Synonyms:

nuclear factor of κ ; light polypeptide gene enhancer in B-cells inhibitor, alpha; NFKBIA; IKBA; MAD-3; NFKBI ;

Immunogen:

Human NFKBIA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20 $^{\circ}$ C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A