

# NFKBIB Antibody



PACO10837

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

NFKBIB

**Uniprot**

Q15653

**Synonyms:**

nuclear factor of  $\kappa$  light polypeptide gene enhancer in B-cells inhibitor, beta; NFKBIB; IKBB; TRIP9 ;

**Immunogen:**

Human NFKBIB.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20 $^{\circ}$ C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A