

NFS1 Antibody



PACO10838

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NFS1

Uniprot

Q9Y697

Synonyms:

NFS1 nitrogen fixation 1 homolog (S. cerevisiae); NFS1; HUSSY-08; IscS; NIFS ;

Immunogen:

Human NFS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A