

NFX1 Antibody



PACO10839

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NFX1

Uniprot

Q12986

Synonyms:

nuclear transcription factor, X-box binding 1; NFX1; DKFZp779G2416; MGC20369; NFX2 ;

Immunogen:

Human NFX1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A