

NGDN Antibody



PACO10843

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NGDN

Uniprot

Q8NEJ9

Synonyms:

neuroguidin, EIF4E binding protein; NGDN; C14orf120; DKFZp564O092; LCP5; NGD; lpd-2 ;

Immunogen:

Human NGDN.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A