

PACO10845

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NGF

**Uniprot**

P01138

**Synonyms:**

nerve growth factor (beta polypeptide); NGF; Beta-NGF; HSAN5; MGC161426; MGC161428; NGFB ; proNGF

**Immunogen:**

Human NGF.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A