

NGFR Antibody



PACO10846

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NGFR

Uniprot

P08138

Synonyms:

nerve growth factor receptor (TNFR superfamily, member 16); NGFR; CD271; Gp80-LNGFR; TNFRSF16; p75(NTR); p75NTR ;

Immunogen:

Human NGFR.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A