

# NHLRC1 Antibody



PACO10849

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NHLRC1

**Uniprot**

Q6VVB1

**Synonyms:**

NHL repeat containing 1; NHLRC1; EPM2A; EPM2B; MALIN; MGC119262; MGC119264; MGC119265; bA204B7.2 ;

**Immunogen:**

Human NHLRC1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A