

NIF3L1 Antibody



PACO10852

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NIF3L1

Uniprot

Q9GZT8

Synonyms:

NIF3 NGG1 interacting factor 3-like 1 (S. pombe); NIF3L1; ALS2CR1; CALS-7; MDS015 ;

Immunogen:

Human NIF3L1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A