

PACO10855

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NIPBL

**Uniprot**

Q6KC79

**Synonyms:**

Nipped-B homolog (Drosophila); NIPBL; CDLS; CDLS1; DKFZp434L1319; FLJ11203; FLJ12597; FLJ13354; FLJ13648; FLJ44854; IDN3; IDN3-B ;

**Immunogen:**

Human NIPBL

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A