

NIPSNAP3A Antibody



PACO10856

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

NIPSNAP3A

Uniprot

Q9UFN0

Synonyms:

nipsnap homolog 3A (C. elegans); NIPSNAP3A; DKFZp564D177; FLJ13953; HSPC299; MGC14553; NIPSNAP4; TASSC ;

Immunogen:

Human NIPSNAP3A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A