

PACO10864

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

NKTR

**Uniprot**

P30414

**Synonyms:**

natural killer-tumor recognition sequence; NKTR; DKFZp686F1754; DKFZp686G0426; DKFZp686J06106; DKFZp686N24126; MGC90527; p104 ;

**Immunogen:**

Human NKTR.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A