

# NKX2-5 Antibody



PACO10865

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NKX2-5

**Uniprot**

P52952

**Synonyms:**

NK2 transcription factor related, locus 5 (Drosophila); NKX2-5; CHNG5; CSX; CSX1; NKX2.5; NKX2E; NKX4-1 ;

**Immunogen:**

Human NKX2-5.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A