

NKX3-1 Antibody



PACO10866

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NKX3-1

Uniprot

Q99801

Synonyms:

234aa,26kd;

Immunogen:

Human NKX3-1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A