

# NLRX1 Antibody



PACO10871

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

NLRX1

**Uniprot**

Q86UT6

**Synonyms:**

NLR family member X1; NLRX1; CLR11.3; DLNB26; FLJ21478; MGC131937; MGC21025; NOD26; NOD5; NOD9 ;

**Immunogen:**

Human NLRX1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A