

NMD3 Antibody



PACO10872

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NMD3

Uniprot

Q96D46

Synonyms:

NMD3 homolog (S. cerevisiae); NMD3; CGI-07; FLJ21053 ;

Immunogen:

Human NMD3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A