

NME2 Antibody



PACO10874

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NME2

Uniprot

P22392

Synonyms:

non-metastatic cells 2, protein (NM23B) expressed in; NME2; MGC111212; NDPK-B; NDPKB; NM23-H2; NM23B; puf ;

Immunogen:

Human NME2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A