

NOL12 Antibody



PACO10887

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

NOL12

Uniprot

Q9UGY1

Synonyms:

nucleolar protein 12; NOL12; FLJ34609; MGC3731; Nop25; dJ37E16.7 ;

Immunogen:

Human NOL12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A