

NOP56 Antibody



PACO10895

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NOP56

Uniprot

O00567

Synonyms:

NOP56 ribonucleoprotein homolog (yeast); NOP56; NOL5A ;

Immunogen:

Human NOP56.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A