

NOTCH3 Antibody



PACO10900

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NOTCH3

Uniprot

Q9UM47

Synonyms:

Notch homolog 3 (Drosophila); NOTCH3; CADASIL; CASIL ;

Immunogen:

Human NOTCH3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A