

NOVA2 Antibody



PACO10901

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NOVA2

Uniprot

Q9UNW9

Synonyms:

neuro-oncological ventral antigen 2; NOVA2; ANOVA; NOVA3 ;

Immunogen:

Human NOVA2.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C,
Avoid freeze / thaw cycles.

Product Images

N/A

N/A