

NOX4 Antibody



PACO10902

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NOX4

Uniprot

Q9NPH5

Synonyms:

NADPH oxidase 4; NOX4; KOX; KOX-1; RENOX ;

Immunogen:

Human NOX4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A