

# NPTX1 Antibody



PACO10912

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NPTX1

**Uniprot**

Q15818

**Synonyms:**

neuronal pentraxin I; NPTX1; DKFZp686J2446; MGC105123; NP1 ;

**Immunogen:**

Human NPTX1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A