

# NR4A2 Antibody



PACO10925

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NR4A2

**Uniprot**

P43354

**Synonyms:**

nuclear receptor subfamily 4, group A, member 2; NR4A2; HZF-3; NOT; NURR1; RNR1; TINUR ;

**Immunogen:**

Human NR4A2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A