

NR5A1 Antibody



PACO10926

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NR5A1

Uniprot

Q13285

Synonyms:

nuclear receptor subfamily 5, group A, member 1; NR5A1; AD4BP; ELP; FTZ1; FTZF1; SF-1; SF1 ;

Immunogen:

Human NR5A1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A