

NRXN1 Antibody



PACO10937

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NRXN1

Uniprot

Q9ULB1

Synonyms:

neurexin 1; NRXN1; DKFZp313P2036; FLJ35941; Hs.22998; KIAA0578 ;

Immunogen:

Human NRXN1.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A