

NTRK3 Antibody



PACO10946

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NTRK3

Uniprot

Q16288

Synonyms:

neurotrophic tyrosine kinase, receptor, type 3; NTRK3; TRKC; gp145(trkC) ;

Immunogen:

Human NTRK3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A