

NUAK2 Antibody



PACO10947

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NUAK2

Uniprot

Q9H093

Synonyms:

NUAK family, SNF1-like kinase, 2; NUAK2; DKFZp434J037; DKFZp686F01113; FLJ90349; SNARK ;

Immunogen:

Human NUAK2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A