

NUDCD2 Antibody



PACO10953

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

NUDCD2

Uniprot

Q8WVJ2

Synonyms:

NudC domain containing 2; NUDCD2; DKFZp686E10109 ;

Immunogen:

Human NUDCD2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A