

NUDT9 Antibody



PACO10961

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NUDT9

Uniprot

Q9BW91

Synonyms:

nudix (nucleoside diphosphate linked moiety X)-type motif 9; NUDT9; MGC3037; NUDT10 ;

Immunogen:

Human NUDT9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A