

NUMB Antibody



PACO10965

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NUMB

Uniprot

P49757

Synonyms:

numb homolog (Drosophila); NUMB; S171 ;

Immunogen:

Human NUMB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A