

NUSAP1 Antibody



PACO10979

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NUSAP1

Uniprot

Q9BXS6

Synonyms:

nucleolar and spindle associated protein 1; NUSAP1; ANKT; BM037; FLJ13421; LNP; PRO0310p1; Q0310; SAPL ;

Immunogen:

Human NUSAP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A