

OGG1 Antibody



PACO10999

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

OGG1

Uniprot

O15527

Synonyms:

8-oxoguanine DNA glycosylase; OGG1; HMMH; HOGG1; MUTM; OGH1 ;

Immunogen:

Human OGG1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A