

OXER1 Antibody



PACO11034

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

OXER1

Uniprot

Q8TDS5

Synonyms:

oxoeicosanoid (OXE) receptor 1; OXER1; GPCR; GPR170; TG1019 ;

Immunogen:

Human OXER1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A