

P2RX4 Antibody



PACO11037

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

P2RX4

Uniprot

Q99571

Synonyms:

purinergic receptor P2X, ligand-gated ion channel, 4; P2RX4; P2X4; P2X4R ;

Immunogen:

Human P2RX4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A