

P2RX7 Antibody



PACO11039

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

P2RX7

Uniprot

Q99572

Synonyms:

purinergic receptor P2X, ligand-gated ion channel, 7; P2RX7; MGC20089; P2X7 ;

Immunogen:

Human P2RX7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A