

P4HA1 Antibody



PACO11042

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

P4HA1

Uniprot

P13674

Synonyms:

prolyl 4-hydroxylase, alpha polypeptide I; P4HA1; P4HA ;

Immunogen:

Human P4HA1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A