

P4HA2 Antibody



PACO11043

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

P4HA2

Uniprot

O15460

Synonyms:

prolyl 4-hydroxylase, alpha polypeptide II;

Immunogen:

Human P4HA2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A