

# PA2G4 Antibody



PACO11044

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

PA2G4

**Uniprot**

Q9UQ80

**Synonyms:**

proliferation-associated 2G4, 38kDa; PA2G4; EBP1; HG4-1; p38-2G4 ;

**Immunogen:**

Human PA2G4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A