

PABPC3 Antibody



PACO11045

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PABPC3

Uniprot

Q9H361

Synonyms:

poly(A) binding protein, cytoplasmic 3; PABPC3; PABP3; PABPL3 ;

Immunogen:

Human PABPC3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A