

# PAFAH1B3 Antibody



PACO11052

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PAFAH1B3

**Uniprot**

Q15102

**Synonyms:**

26kd;

**Immunogen:**

Human PAFAH1B3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A