

PACO11054

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PAICS

Uniprot

P22234

Synonyms:

phosphoribosylaminoimidazole carboxylase, phosphoribosylaminoimidazole succinocarboxamide synthetase; PAICS; ADE2; ADE2H1; AIRC; DKFZp781N1372; MGC1343; MGC5024; PAIS ;

Immunogen:

Human PAICS.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A