

PAK4 Antibody



PACO11059

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

PAK4

Uniprot

O96013

Synonyms:

p21 protein (Cdc42/Rac)-activated kinase 4;

Immunogen:

Human PAK4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A