

PAOX Antibody



PACO11067

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PAOX

Uniprot

Q6QHF9

Synonyms:

polyamine oxidase (exo-N4-amino); PAOX; DKFZp434J245; MGC45464; PAO ;

Immunogen:

Human PAOX.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A