

PARP1 Antibody



PACO11073

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PARP1

Uniprot

P09874

Synonyms:

poly (ADP-ribose) polymerase 1; PARP1; ADPRT; ADPRT1; PARP; PARP-1; PPOL; pADPRT-1 ;

Immunogen:

Human PARP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A