

PARP11 Antibody



PACO11074

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PARP11

Uniprot

Q9NR21

Synonyms:

poly (ADP-ribose) polymerase family, member 11; PARP11; C12orf6; DKFZp779H0122 ;

Immunogen:

Human PARP11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A